




**DENTAL
INSTRUMENTS
CATALOGUE**

 **DENTOGRADE**
MANUFACTURERS

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision there Grooves for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 20 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours.



- High-polish mirror finish tips.
- Manufactured within the guidelines of the FDA quality assurance system for good manufacturing practices.
- Highest grade German stainless steel, razor sharp, durable tips. Rockwell hardened to 60 RHc.
- 27 of the most popular designs from graceys to McCalls to columbias, and many european patterns.
- All instruments are guaranteed against defects in material and workmanship.

- All instruments bear the ce mark, emblematic of the european community's highest standards.
- Customer satisfaction guarantee, too!
- Less than 20 grams ... Lightweight for less hand fatigue.
- Manufactured to the most exacting standards so the 50th instrument you buy is exactly like your first.
- You decide ... Why pay more for the same high-quality instruments?

Dear Valued Buyers

We feel immense pleasure to present our latest catalogue in which we have featured a most comprehensive range of all types of quality Surgical & Dental instruments required for day to day use for Hospital. We will be pleased to quote you a competitive Price and Delivery time.

All instruments are manufactured according to the recognized standard patterns currently in use in Hospital throughout the world and are being produced under restrict control of experts to avoid any flaw and to refrain any complaint at any stage.

A group of very vast experienced persons in manufacturing handling our quality Control System.

Any instruments proven defective in workmanship or material will be replaced with out charges.

We express our thanks to all our customers for the confidence they have placed us in the past. We further assure them that we will continue our best efforts to gain their confidence in the future.

We are always at your service.



Chief Executive

**WE BRING TECHNICAL KNOWLEDGE & EXPERIENCE AT ONE PLACE.
WE INCORPORATE TECHNOLOGY FOR BETTER PRODUCTIVITY.**



CNC LATHE WORKING MACHINE

Hardness Test
Machine

CNC Milling MACHINE WORKING

MANUFACTURING PROCESS

HANDLES FOR MOUTH MIRRORS, TWEEZERS

Ручки для зеркал
Manches pour miroirs à bouche
Mundspiegelgriffe
防盗镜手柄

Ручки для дзеркал
Mangos para espejos de boca
مقابض لمرآيا الفم




PLANE/RHODIUM FRONT SURFACE

DOUBLE SIDE MOUTH MIRRORS RHODIUM



DG1-1036 Mirror Handle Solid Stem with Rule 35mm



Зонды
Sondes
Zahnsonden
探索者

Зонди
Sondas
المستكشفين

DG1-1047 Explorer # 23 SE

DG1-1048 Explorer # 6 SE

DG1-1049 Explorer # 8A SE

DG1-1050 Explorer #5 DE

DG1-1051 Explorer DG 16 DE

DG1-1052 Explorer 23/6 DE

DG1-1053 CP-C

DG1-1054 CP-11


DG1-1055 CP-12 without Ball

DG1-1056 CP-12 with Ball Marking

DG1-1057 CP-12S O.M.S.


DG1-1058 Williams

DG1-1059 Goldman Fox



Зонды
Sondes
Zahnsonden
探索者

Зонди
Sondas
المستكشفين






Available in Hollow Handle 8mm, Ø9.5mm

DG1-1082 S Gracey Curette Set of 9 Pcs

Gracey Curette Gold Coated Pointed Set of 7 Pcs




Gracey Curette Standard
Originally designed to be finishing curettes for fine scaling and root planing.

Each Gracey blade is offset at 70° and has one cutting edge - the lower edge. The combination of this unique blade and 9 different shank designs for specific teeth surface provide improved adaptation and deposit removal in scaling procedures.



DG1-1092
Gracey Curette Gold Coating)
Set of 7 Pcs



DG1-1100 Set of 6 Pcs


GRACEY MINI FIVE

Available in Hollow Handle 8mm, Ø9.5mm




DG1-1102 Set of 6 Pcs

GRACEY AFTER FIVE

Available in Hollow Handle 8mm, Ø9.5mm




Gracey Curette Mini
Terminal shank is 3mm longer than Gracey Standard compared to the standard Gracey curette, blade is 50% shorter and 10% thinner
Good for scaling in deep and/or narrow pockets



CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس
 水果

SICKLE SCALERS
Available in Hollow Handle 8mm, Ø9.5mm




Кюрети
Curetas
Küretten
水果


Кюрети
Curette
کرتس

Available in Hollow Handle Ø8.0mm - 9.5mm

DG1-1134 SET
Sugarman File Set



Spoon Curettes Oval Set of 3 Pcs




DG1-1146 SET OF 10 PCS




DG1-1147
Ochsenbein
Periodontal Chisel



DG1-1148
Ochsenbein
Periodontal Chisel



Williger Curette
13.5 cm




QUALITY INNOVATION

EXCAVATORS

Экскаваторы
Excavateurs
Excavatoren

Экскаватори
Excavadores
حفارات

Available in Hollow Handle 8mm, Ø9.5mm



■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ 填充仪器

■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك

Available in Hollow Handle Ø8mm - 9.5mm

DG1-1178 SET COMPOSITE GOLD TIP INSTRUMENTS SET OF 8 PCS



DG1-1178 Double Paddle #1T 7 mm long/2.0 mm wide.



DG1-1179 Paddle/Condenser #3T 7mm long/2mm wide with a 1.1mm



DG1-1180 Double Paddle
#4T 11mm long/2.0mm wide.



DG1-1183 Paddle/Acorn
#23T 11mm long/2mm wide 2.3mm diameter.



DG1-1181 Paddle/Condenser
#20T 11mm long/2mm wide with a 2.0mm



DG1-1184 Duck-Head
#28T 3.8mm 2.5 mm.



DG1-1182
Carver Hollenback 3



DG1-1185
Carver Hollenback 3S



DG1-1186 SET PLUGGER, CONDENSER & BURNISHER SET OF 6 PCS



DG1-1186 Multifunction Plugger & Condenser PLGS/M Plugger Medium



DG1-1187
Multifunction Plugger & Condenser
PLGS/L Plugger Large



DG1-1189
Universal Amalgam Kto Titanium Coated



DG1-1188 B27/29 BALL
BURNISHER - Titanium Coated



DG1-1190
EGG Burnisher DE



DG1-1191 T Ball Burnisher

FILLING INSTRUMENTS


■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ 填充仪器

■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك

Available in Hollow Handle 8mm, Ø9.5mm


DG1-1192 SET

COMPOSITE GOLD TIP INSTRUMENTS SET OF 7 PCS




DG1-1199 SET

COMPOSITE BLACK TITANIUM COLOUR COATED INSTRUMENTS SET OF 5 PCS



DG1-1204 SET MARGIN TRIMMERS SET OF 6 PCS



FILLING INSTRUMENTS

■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ 填充仪器

■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك

DG1-1210

P.K. THOMAS WAX CARVER SET OF 5 PCS



Available in Hollow Handle Ø8.0mm - 9.5mm

COMPOSITE INSTRUMENTS BLUE COATED



DG1-1212

Composite Instruments
Set of 6 Pcs

Blue Color Coated with Cassettes



DG1-1211

P.K. THOMAS WAX CARVER
STAINLESS STEEL
SET OF 5 PCS



DG1-1213 A SET

Composite Filling
Instruments Blue Coated
Without Handle Color
Set of 6 Pcs

DG1-1214

Composite Filling
Instruments Gold Coated
Set of 6 Pcs

FILLING INSTRUMENTS

■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ 填充仪器

■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك



Available in Hollow Handle Ø8.0mm - 9.5mm

COMPOSITE INSTRUMENTS SET OF 5 PCS


DG1-1215 SET



Available in Hollow Handle Ø8.0mm - 9.5mm

COMPOSITE INSTRUMENTS SET OF 5 PCS


DG1-1216 SET B



Silicon Handle

COMPOSITE INSTRUMENTS SET OF 5 PCS

DG1-1217



FILLING INSTRUMENTS

■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ 填充仪器

■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك

Silicon Handle

Eccesso

Misura

Replica Posterior

Posterior Misura

Cusp Misura

Fissura

Modella

Condensa

Silicon Handle

Posterior Misura

Cusp Misura

Fissura

Modella


Condensa



DG1-1218
Set of 8 Pcs




DG1-1219
Composite
Set of 5 Pcs



■ # 1/2 - # 3/4 - # 5/6
 ■ # 7/8 - # 9/10
 ■ # 11/12 - # 15/16
 ■ # 13/14 - # 17/18



DG1-1220
Composite
Set of 7 Pcs



■ # 1/2 - # 3/4 - # 5/6
 ■ # 7/8 - # 9/10
 ■ # 11/12 - # 15/16
 ■ # 13/14 - # 17/18

Ножи для воска
Couteaux a cire
Wachsmesser
蜡刀

Ножі для воску
Cuchillos para cera
السكاكين الشمع

WAX KNIVES



DG1-1221
Fig. 1
Fahnenstock
125mm



DG1-1222
Fig. 2
Fahnenstock
170mm



DG1-1223
Fig. 1
Lessmann
125mm



DG1-1224
Fig. 2
Lessmann
170mm



DG1-1225
Fig. 1
with small spoon
125mm



DG1-1226
Fig. 2
with large spoon
170mm



DG1-1227
Plaster Knife

WAX SPATULAS



DG1-1228
Wooden Handle
19.5cm curved



DG1-1229
Wooden Handle
20cm curved



DG1-1230
Wooden Handle
22.5cm curved



DG1-1231
Wooden Handle
19.5cm



DG1-1232
Wooden Handle
21.5cm



DG1-1233
Wooden Handle
20.5cm



DG1-1234
Fig. 3
12.5cm



DG1-1235
Fig. 3
12.5cm



DG1-1236
Fig. 4
13.5cm

DG1-1237
Carvers Set of 10 Pcs



DG1-1238
Fig. 1



DG1-1239
Fig. 2



DG1-1240
Fig. 3



DG1-1241
Graded Handle
Hylin Carver



DG1-1242
Hylin



DG1-1243
Hylin-Le Cron



DG1-1244
Hylin-Ceramica



DG1-1245
Hylin



DG1-1246
Le Cron



DG1-1247
Le Cron



DG1-1248
Le Cron



DG1-1249
Le Cron

Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (نمط اللغة الإنجليزية)

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

DG1-1250 SET

LUXATOR SET OF 4 PCS STRAIGHT



DG1-1250

DG1-1251

DG1-1252

DG1-1253



DG1-1254 SET

LUXATOR SET OF 4 PCS CURVED



DG1-1254

DG1-1255

DG1-1256

DG1-1257

APICAL ROOT ELEVATORS



DG1-1258
Fig. 1, 3mm
Frideman
Straight



DG1-1259
Fig. 2
Frideman
Left



DG1-1260
Fig. 3
Frideman
Right



lower roots elevator



DG1-1261
Fig. 304



DG1-1262
Apical 302A



DG1-1263
Apical 303A

DG1-1263 A ROOT ELEVATOR SET OF 10 PCS



DG1-1264
Root Elevator
DL 32
for upper roots
round point.



DG1-1265
Root Elevator
DL 33
for upper roots
round point.



DG1-1266
Bein
DL 34
for upper roots
round point.



DG1-1267
Flohr
Fig. 3
for lower roots



DG1-1268
Flohr
Fig. 4
for lower roots



DG1-1269
Cryer-White
for lower roots,
wide pattern



DG1-1270
Cryer-White
for lower roots,
wide pattern



DG1-1271
Apical 301



DG1-1272
Apical 302



DG1-1273
Apical 303

Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (نمط اللغة الإنجليزية)

This section includes the Luxating Elevator set with Stainless Steel Handle.
Long & light weight handle design provide excellent smooth grip with controlled force transfer.
The working procedure is similar as mentioned in the previous section.

Luxating Elevators Set



DG1-1280A BEIN ROOT ELEVATORS SET OF 6 PCS



DG1-1286A BEIN ROOT ELEVATORS SET OF 6 PCS



ROOT ELEVATORS

 Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)
 提取强制 (英语模式)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج (نمط اللغة الإنجليزية)

Wisdom Roots Elevator




DG1-1293
Wisdom
for 18 and 38

DG1-1294
Wisdom
for 28 and 48

Potts Elevators

These root elevators are used for releasing the ligaments around the neck of the tooth, and to facilitate the extraction. The working end are smooth. These are used for impacted upper 3rd Molars without bony obstruction.

Elevator according to Potts
Working Tips are smooth
Smooth Surface Handle
Surgical Stainless Steel
Sterilized in all Commercial
Sterilizers



6 Potts Elevator

7 Potts Elevator

DG1-1295

DG1-1296

Schmeckebier Apexo Elevators

Schmeckebier Apexo elevators use to loosen the tooth from the periodontal ligament and ease extraction. Sand Blasted no slip handle. Surgical Stainless Steel. Sterilized in all Commercial Sterilizers.



Schmeckebier Apexo Elevator L

Schmeckebier Apexo Elevator R

DG1-1297

DG1-1298

Apical Elevators

Apical root elevators, right, left. For upper and lower root tips. Width of tip 3 mm. Smooth Surface Handle. Surgical Stainless Steel. Sterilized in all Commercial Sterilizers.



Apical elevator Left, 3mm

Apical elevator Right, 3mm

DG1-1299

DG1-1300

Seldin Elevators

Seldin root elevators. Smooth Surface Handle. Surgical Stainless Steel. Sterilized in all Commercial Sterilizers.



34 Str. Seldin Elevator, Upper and lower anterior roots

34S Seldin Elevator, Upper & lower anterior roots, sharp cone edge

25 Seldin Elevator, Lower 3rd Molars

DG1-1301

DG1-1302

DG1-1303

B-Cogswell & Apexo Elevators

Smooth Surface Handle. Surgical Stainless Steel. Sterilized in all Commercial Sterilizers.



301 Apexo Elevator, for displacement of deeply seated roots


B Cogswell Elevator

DG1-1304

DG1-1305

Medium Cryer Root Elevators set

Left & Right, with medium working points. Working Tips are smooth. Smooth Surface Handle. Surgical Stainless Steel. Sterilized in all Commercial Sterilizers.




21Mini Cryer Elevator, Right

22Mini Cryer Elevator, Left

DG1-1306

DG1-1307

Left & Right, with large working points. Working Tips are smooth. Smooth Surface Handle. Surgical Stainless Steel. Sterilized in all Commercial Sterilizers.




44 Cryer Root Elevator, Right

44 Cryer Root Elevator, Left

DG1-1308

DG1-1309



30 Cryer Root Elevator, Right

31 Cryer Root Elevator, Left

DG1-1310

DG1-1311

ROOT ELEVATORS

Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج (نمط اللغة الإنجليزية)

Posterior Flared Tip Elevators

This section includes a comprehensive set of Posterior Flared Tip Elevators. These elevators used to loosen the tooth from the periodontal ligament and ease extraction, particularly used for Upper Posterior Roots and Upper 3rd Molar. E46 & E46R: Mesial & Distal with 2.75mm, Smooth working tips E77 & E77R: Mesial & Distal with 3.75mm, Serrated working tips Surgical Stainless Steel Hollow handle, very light weight Completely Sterilizable



DG1-1312
46R Elevator, Smooth,
Mesial, 2.75mm



DG1-1313
46L Elevator, Smooth,
Distal, 2.75mm



DG1-1314
77R Elevator, Serrated,
Mesial, 3.75mm



DG1-1315
77 Elevator, Serrated,
Distal, 3.75mm

Lindo-Levien Elevators

These root elevators are used for releasing the ligaments around the neck of the tooth, and to facilitate the extraction. The working end serration can be very helpful during the initial phase of extraction.

Elevator according to Lindo-Levien
Working part straight with serration
Width 2.2mm 3.0mm or 4.0mm
Smooth Surface Handle
Surgical Stainless Steel
Sterilized in all Commercial Sterilizers



Upper root elevator with helicoidal tip

lower roots elevator



DG1-1316
Lindo-Levien
Elevators
3mm, Straight



DG1-1317
Lindo-Levien
Elevators
3mm, Curved




DG1-1318
Lindo-Levien
Elevators
4mm, Straight



DG1-1319
Lindo-Levien
Elevators
4mm, Curved

LENDO-LEVIAN ELEVATORS STR.




DG1-1320

DG1-1321

DG1-1322



LENDO-LEVIAN ELEVATORS CVD.



DG1-1323
Fig. 1

DG1-1324
Fig. 2

DG1-1325
Fig. 3



ROOT ELEVATORS


Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (تمط اللغة الإنجليزية)


DG1-1326 LUXATORS SET OF 7 PCS



DG1-1334 LUXATORS SET OF 10 PCS



BEIN ROOT ELEVATORS SET OF 6 PCS



DG1-1351
Bein Root Elevators
Set of 6 Pcs

ROOT ELEVATORS

Шипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (نمط اللغة الإنجليزية)




DG1-1359
Root Elevator Set of 7 Pcs

Luxator Root Elevators




FLOHR ROOT ELEVATORS SET OF 6 PCS



ROOT ELEVATORS


Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (تمط اللغة الإنجليزية)



DG1-1385
Root Elevator Set of 7 Pcs

ROOT ELEVATORS SET OF 7 PCS



DG1-1393
Root Elevator Set of 7 Pcs
with Gold handle



DG1-1399
Root Elevators
Set of 5 Pcs



DG1-1404
Root Elevators
Set of 4 Pcs

ROOT ELEVATORS

Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (نمط اللغة الإنجليزية)



Root Pick Elevator With Gold Flexible Points



DG1-1405



DG1-1406



DG1-1407



DG1-1408



DG1-1409



DG1-1410



DG1-1411

Atraumatic Extraction Bayonet Root Elevators



DG1-1412
Bayonet Elevator
Straight



DG1-1413
Bayonet Elevator
Offset "In" Small



DG1-1414
Bayonet Elevator
Offset "Out" Small



DG1-1415
Bernard Bayonet
Elevator Offset Left



DG1-1416
Bernard Bayonet
Elevator Offset Right



DG1-1417
PDL Elevators Set of 7 Pcs



DG1-1418
PDL Elevators Set of 7 Pcs Color Coated

ROOT ELEVATORS

Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (نمط اللغة الإنجليزية)

Heidbrink Roots Elevator



DG1-1419 DG1-1420 DG1-1421

Heidbrink Root tip Picks



DG1-1422
Fig. 1

Syndesmotomes Roots Elevator



DG1-1425 DG1-1426 DG1-1427
Straight Left Right



DG1-1428 DG1-1429 DG1-1430 DG1-1431
Hylin 3mm Str. Overall length 154mm
Hylin 4mm Str.
Hylin 3.5mm Left
Hylin 3.5mm Right



Coupland



universal mini-elevator for fragments



syndesmotome

DG1-1435
Fig. 1
Syndesmotome

HEIDBRINK ROOT TIP SET OF 3 PCS

Available in Hollow Handle 8mm, Ø9.5mm



RETRACTORS - MOUTH GAGS


- | | |
|----------------------------------|--------------------------|
| Роторасширители | Роторозширювачі |
| Ecarteurs - Ouvre-bouches | Separadores - Abre-bocas |
| Vestibulumabhalter / Mundspanner | الكامشات - الكمامات الفم |
| 经销商-嘴塞 | |



DG1-1442
Stainless Steel Single
Ended Cheek Retractors



DG1-1443
Stainless Steel Double
Ended Cheek Retractor



DG1-1444
Cheek Retractor
Stainless Steel Long



DG1-1445
Dental Lip
Retractor



DG1-1446
Left



DG1-1447
Straight
Retractor Blue Titanium Coating



DG1-1448
Right



DG1-1449
Retractor
Blue Titanium Coating



DG1-1450
Stainless Steel Double
Black & Gold Ended Cheek
Retractors

CRURBAN, University of Minnesota modified cheek retractor, is used to hold mucoperiosteal flaps, cheeks, lips and tongue away from the surgical area. 20.6mm



DG1-1451
Urban Cheek
Retractor

The Minnesota Retractor is used to hold mucoperiosteal flaps, cheeks, lips and tongue away from the surgical area. Retractor is 5.5" and 20.5mm wide and is made from surgical grade stainless steel for longevity and corrosion resistance.



DG1-1452
University of Minnesota
Retractor

The Camwood-Minnesota Cheek Retractor is used to hold mucoperiosteal flaps, cheeks, lips and tongue away from the surgical area. 5.9" long. Same design as the CRM Minnesota Retractor with longer handle end. Made with surgical grade stainless steel.



DG1-1453
Camwood Minnesota
Cheek Retractor



DG1-1454
Minnesota



DG1-1455
Minnesota Cheek
Retractor 15.5cm



DG1-1456
Brunings
Cheek Retractor
19cm



DG1-1457
Cheek Retractor
20cm
Easy to Bend
by Hand

RETRACTORS - MOUTH GAGS

Роторасширители
 Ecarteurs - Ouvre-bouches
 Vestibulumabhalter / Mundspanner
 经销商-嘴塞

Роторозширювачі
 Separadores - Abre-bocas
 الكامشات - الكامشات الفم



DG1-1458 10.5cm
DG1-1459 13.5cm
 Weitlaner
 Fig. 1



DG1-1460 10.5cm
DG1-1461 13.5cm
 Weitlaner
 Fig. 2



DG1-1462
 Mouth Gag Molt
 9cm for (Children)



DG1-1463
 Mouth Gag Molt
 11.5cm for (Adults)




DG1-1464
 Mouth Gag Molt
 14cm for (Adults)



DG1-1465
 Cheek Retractor
 Long



DG1-1466
 Cheek Retractor
 Long



DG1-1467
 14cm

DG1-1468
 15cm

Cheek Retractor



DG1-1469
 Unique Retractor
 Length 17.5cm
 angled 45°




DG1-1470
 Unique Retractor
 Length 17.5cm
 angled 90°



DG1-1471
 Cheek Retractor
 Combination
 21cm, 2 Pcs



DG1-1472
 Sternberg
 14 cm
 Double Ended



DG1-1473
 Sternberg
 16 cm



DG1-1474
 Sternberg
 14 cm
 Single Ended

Инструменты для системы коффердам
Pincers porte-crampons a digue
Kofferdam Klammerzangen
橡胶坝夹钳

Инструменты для системы коффердам
Pinzas de engrapar para dique de caucho
السداد المطاط الملقط المشابك



Rubber Dam Kits



DG1-1475

On request the Content of the rubber dam kit can be changed with other instruments.

Content:


- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Rubber Dam Frame L
- 05 Rubber Dam Frame M
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

- 08 Rubber Dam Clamp Fig. 212
- 09 Rubber Dam Clamp Fig. 211
- 10 Rubber Dam Clamp Fig. 8
- 11 Rubber Dam Clamp Fig. 8A
- 12 Rubber Dam Clamp Fig. 14
- 13 Rubber Dam Clamp Fig. 14A
- 14 Rubber Dam Clamp Fig. 0
- 15 Rubber Dam Clamp Fig. 00
- 16 Rubber Dam Clamp Fig. 2
- 17 Rubber Dam Clamp Fig. 2A
- 18 Rubber Dam Clamp Fig. 7



Rubber Dam Kit DG1-1476

Anisworth 17cm, Ivory 17cm, Brewer 17.5cm
Rubber Dam Frame (Children), Rubber Dam Frame (Adult)
Rubber Dam Clamps 11 Pcs
Fig. 00, Fig.0, Fig.1, Fig.2A, Fig.3, Fig.7,
Fig.8, Fig.8A, Fig.9, Fig.212, Fig.14




DG1-1477
Rubber Dam Punch



DG1-1478
Brewer



DG1-1479
Plier For Clamp



DG1-1480 SET
Enc Clamp Kit 11 Pcs
Rubber Dam Clamps 9 Pcs
Fig.210, Fig.211, Fig.2, Fig.2A, Fig.207
Fig.201, Fig.202, Fig.203, Fig.204

RUBBER DAM CLAMPS

Резиновая дамба зажимы комплекты
Trousses de caoutchouc
Gummidammklammern
橡胶夹

Гумова дамба затискачі комплекти
Conjuntos de abrazaderas de retención de caucho
المطاط السد المشابك مجموعات



DG1-1481
Fig. 15



DG1-1482
Fig. 18



DG1-1483
Fig. 24



DG1-1484
Fig. 25



DG1-1485
Fig. 26



DG1-1486
Fig. 28



DG1-1487
Fig. 30



DG1-1488
Fig. 31



DG1-1489
Fig. 51



DG1-1490
Fig. 138



DG1-1491
Fig. 139



DG1-1492
Fig. 200



DG1-1493
Fig. 201



DG1-1494
Fig. 202



DG1-1495
Fig. 203



DG1-1496
Fig. 204



DG1-1497
Fig. 205



DG1-1498
Fig. 3



DG1-1499
Fig. W3



DG1-1500
Fig. 4



DG1-1501
Fig. 56A



DG1-1502
Fig. 7



DG1-1503
Fig. 7A



DG1-1504
Fig. 8



DG1-1505
Fig. 8A



DG1-1506
Fig. W8A



DG1-1507
Fig. 10



DG1-1508
Fig. 11



DG1-1509
Fig. 12A



DG1-1510
Fig. 13A



DG1-1511
Fig. 14



DG1-1512
Fig. 14A



DG1-1513
Fig. 14AN



DG1-1514
Fig. 22



DG1-1515
Fig. 27



DG1-1516
Fig. 207



DG1-1517



DG1-1518
Fig. 206



DG1-1519
Fig. 208



DG1-1520
Fig. 209



DG1-1521
Fig. 1



DG1-1522
Fig. 1A



DG1-1523
Fig. 2



DG1-1524
Fig. 2A



DG1-1525
Fig. W2A



DG1-1526
Fig. 0



DG1-1527
Fig. 00



DG1-1528
Fig. 210



DG1-1529
Fig. 211



DG1-1530
Fig. 212



DG1-1531
Fig. 6



DG1-1532
Fig. 9



DG1-1533
Fig. N06



DG1-1534
Fig. U67



DG1-1535
Fig. P1



DG1-1536
Fig. P2



DG1-1537
Fig. G-1



DG1-1538
Fig. G-2

RUBBER DAM CLAMPS

Резиновая дамба зажимы комплекты
Trausses de caoutchouc
Gummidammklammern
橡胶夹

Гумова дамба затискачі комплекти
Conjuntos de abrazaderas de retención de caucho
المطاط السد المشابك مجموعات



DG1-1539

Stainless Steel Gold Plated Rubber Dam Clamp Tray with 11 Different Clamps
Fig. 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



DG1-1540

Stainless Steel Rubber Dam Clamp Tray with 12 Different Clamps
Fig. 209, 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



DG1-1541

Stainless Steel Rubber Dam Clamp Tray with 8 Different Clamps
Fig. 1, 26, 29, 51, 60, 202, 205, 210



DG1-1542

Stainless Steel Rubber Dam Clamp Tray
(without clamps)



DG1-1543

Stainless Steel
Amalgam Well
large



DG1-1544

Ivory
without shields



DG1-1545

Tofflemire-Junior
45 mm



DG1-1546

Tofflemire-Universal
60 mm



DG1-1547

Matrix Ivory 8
6 mm



AMALGAM INSTRUMENTS

TIP SIZES FOR DOUBLE-END

DG1-1548 - 1.5mm-2.0mm

DG1-1549 - 2.0mm-2.8mm

DG1-1550 - 2.0mm-3.2mm

Резиновая дамба зажимы комплекты
Traousses de caoutchouc
Gummidammklammern
橡胶夹

Гумова дамба затискачі комплекти
Conjuntos de abrazaderas de retención de caucho
المطاط السد المشابك مجموعات



DG1-1551
Metal tip
22cm
Hammer With Wooden Handle



DG1-1552
Metal tip
22cm
Hammer With Wooden Handle



DG1-1553
1
For plaster



DG1-1554
2
For metal



DG1-1555
Fret saw
25cm



DG1-1556
Small
Saw Blades



DG1-1557
Large
Saw Blades



DG1-1558
Ceri Saw



DG1-1559
Plaster Shears
20cm



Crucible Tongs
DG1-1560 25cm
DG1-1561 32cm



DG1-1562
Ceri Saw
18cm



DG1-1563
Blade
0,12x5 mm
Conf. 12 pcs



DG1-1564
Blade
0,25x5 mm
Conf. 12 pcs



DG1-1565
Blade
0,25x3 mm
Conf. 12 pcs

CROWN REMOVER

Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen
皇冠去除

Инструменты для видалення коронок
Arranca-corona
مزيلات التاج



Fitting is possible on both ends

Fitting is possible on both ends

DG1-1566

Crown Remover,
For different level
of intensity



The inner spring can be remove, allow to control the stroke manually

DG1-1567

Crown Remover,
For different level
of intensity



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Fig. 5



Fig. 6



A

B

DG1-1568

Crown Remover
Set of 3 Pcs



DG1-1569

Crown Remover
Set of 5 Pcs



B

C

D

E

DG1-1570

Crown Remover
Set of 6 Pcs



Turkish

A

B

C

D

DG1-1571

Crown Remover
Set of 4 Pcs



DG1-1572

Offset Crown
Splitter Remover
Elevator



DG1-1573

Fig. 1



DG1-1574



DG1-1575

Crown Spreaders
Fig. 3



DG1-1576

Crown Spreaders



DG1-1577

Crown Removal Pliers
14cm



Rubber jaws
removes
Crown without
Damage
Crown Remover
Replacement
Soft Pads
DG1-1578 TIP


Free 10 Crown Remover Replacement Soft Pads

 Костные кусачки
 Hohlmeißelzangen
 Pincers-gouges

 Кусачки для кісток
 Pinzas-gubias
 روتنجورس العظام

 Шприцы стоматологические карпульные
 Seringues à usage dentaire
 Dentspritzen
 牙科注射器

 Шприці стоматологічні карпульні
 Jeringas dentales
 المحاقن الأسنان



DG1-1579
 Micro-Friedman 15cm
 A forma "S"
 "S" Shaped

DG1-1580
 Micro-Friedman
 16cm

DG1-1581
 Friedman
 14 cm
 30°

DG1-1582
 Friedman
 14 cm
 45°

DG1-1583
 Micro-Friedman
 14cm -90°

DG1-1584
 Ruskin Mini, Rongeurs
 Curved, 15cm,
 5mm diameter tip

DG1-1585
 Boehler, Rongeurs
 double action hinged, curved,
 15cm, 2mm diameter tip

DG1-1586
 Bayer, Rongeurs
 double action hinged,
 18cm, 3mm diameter tip



Type C
DG1-1587 1.8ml
DG1-1588 2.2ml



3-Ring Non Folding
DG1-1596 1.8ml
DG1-1597 2.2ml



DG1-1602 1.8ml
DG1-1603 2.2ml
 Angled



DG1-1589
 Aspirated Dental Cartridge
 Syringe 1.8 ml



Aspirating Syringe
DG1-1598 1.8mm
DG1-1599 2.2mm



Aluminium - Dark
DG1-1604 1.8 ml
DG1-1605 2.2 ml



DG1-1590 1.8 ml
DG1-1591 2.2 ml
 Aspirated Dental Cartridge
 Syringe Folding Style



Cartridge Syringe Folding
DG1-1600 1.8mm
DG1-1601 2.2mm



Aluminium - Blue
DG1-1606 1.8 ml
DG1-1607 2.2 ml



Projection System
DG1-1592 1.8ml
DG1-1593 2.2ml



Stieglitz Root Splinter Forceps




DG1-1594
 # 1001
 Intraligamental
 Anesthesia
 Syringe in
 Black Colour



DG1-1595 1.8ml
 Pistol-Ject
 injects 0.20ml with every
 click Lever action for
 high injection pressure


Autoclavable



DG1-1608
 Stieglitz
 Root Fragments



DG1-1609
 Stieglitz
 Root Fragments



DG1-1610
 Stieglitz Splinter Forceps
 14cm 90°



DG1-1611
 Splinter Forceps
 Straight
 18cm




DG1-1613
 Splinter Forceps
 18cm 45°

SCISSORS

■ Ножницы десневые
 ■ Ciseaux a gengives
 ■ Zahnfleischscheren
 ■ 口香糖剪刀

■ Ножницы для ясен
 ■ Tijeras para encias
 ■ مقص صمغ




Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren
口香糖剪刀

Ножницы для ясен
Tijeras para encias
مقص صمغ



DG1-1632
Iris Scissors T.C Super Cut
Black Titanium Coated
large rings, curved, 11.5cm



DG1-1633
Goldman-Fox
Scissors,
Standard 13cm



DG1-1634
Goldman-Fox
Scissors,
Super Cut
13cm




DG1-1635
Goldman-Fox
Scissors, T.C
13cm




DG1-1636
Goldman-Fox Scissors,
T.C Super Cut Curved
Black Titanium Coated
13cm



DG1-1637
Goldman-Fox Scissors,
T.C Super Cut Straight
Black Titanium Coated
13cm



DG1-1638
Goldman-Fox
Scissors, Standard
Double curved, 13cm



DG1-1639
Goldman-Fox
Scissors, Super Cut
Double curved, 13cm




DG1-1640
16cm
1 Blade
Serrated
Locklin
Dissecting
Scissors,
S/S




DG1-1641
16cm
1 Blade
Serrated
Locklin
Dissecting
Scissors,
S/S



DG1-1642
La Grange Scissors,
Standard
double curved, 11.5cm



DG1-1643
La Grange Scissors,
Super Cut
double curved, 11.5cm



DG1-1644
La Grange Scissors,
Super Cut
Black Titanium Coated
double curved, 11.5cm



DG1-1645
De Bakey
Dissecting
Scissors
Angled, S/S



DG1-1646
Joseph Scissors,
Standard
curved, 14cm



DG1-1647
Joseph Scissors,
Super Cut
curved, 14cm



DG1-1648
Joseph Scissors,
T.C Super Cut
Black Titanium Coated
curved, 14cm



Potts Dissecting
Scissors
Angled, S/S



25°
DG1-1649
19cm
DG1-1650
23cm



45°
DG1-1651
19cm
DG1-1652
23cm



60°
DG1-1653
18.5cm
DG1-1654
23cm

Иглодержатели
 Porte-aiguilles
 Nadelhalter

Голкотримачі
 Porta-agujas
 أصحاب إبرة

MICRO NEEDLE HOLDER/SCISSORS, FORCEPS AND SCALPEL HANDLES



Forceps
DG1-1655 14cm
DG1-1656 16cm
DG1-1657 18cm



Castroviejo Scissors TC
 Straight
DG1-1658 14cm
DG-1659 18cm




Castroviejo Scissors TC
 Curved
DG1-1660 14cm
DG-1661 18cm



Castroviejo Needle TC
 Straight
DG1-1662 14cm
DG1-1663 18cm



Castroviejo Needle TC
 Curved
DG1-1664 14cm
DG1-1665 18cm



DG1-1666
 Scalpel Handle
 movable



DG1-1667
 Scalpel Handle
 Auto Push Blade



DG1-1668
 Scalpel Handle
 straight, #3



DG1-1669
 Micro Surgical
 Scalpel Handle



DG1-1670
 Surgical Bone
 Files-Miller
 Bf45



DG1-1671
 Surgical Bone
 Files-Miller
 Bf1x



DG1-1672
 Surgical Bone
 Files-Miller
 Bf3x



DG1-1673
 Surgical Bone
 Files-Miller
 Bf21

BONE FILES

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачі
Porta-agujas
اصحاب ابرة

T.C GOLD



EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج المنقط (نمط اللغة الإنجليزية)



Fig. 1
upper laterals
and canines



Fig. 7
upper
premolars



Fig. 13
lower
premolars



Fig. 17
upper molars,
right



Fig. 18
upper molars,
left



Fig. 22
lower roots



Fig. 51A
upper roots



Fig. 67
upper third
molars



Fig. 74N
lower molar



Fig. 79
lower third
molars



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



DG1-1709
Set of 10 Pieces
Fig. 1
Fig. 7
Fig. 13
Fig. 17
Fig. 18
Fig. 22
Fig. 51A
Fig. 67
Fig. 74N
Fig. 79



160mm
Witzel Forceps
Lower roots, Universal
拔除上下颌牙根使用



165mm
Witzel Forceps
Upper roots, Universal
拔除上下颌牙根使用



165mm
Witzel Forceps
Lower roots, Universal
拔除上下颌牙根使用



DG1-1710
Witzel Forceps
With Gold
Set of 3 Pcs

EXTRACTING FORCEPS (For Children)

Щипцы экстракционные детские

Daviers - forme anglaise

Pinzeper per estrazione per bambini

Щипці екстракційні дитячі

Fórceps de extracción - forma inglesa

ستخراج الملقط (للأطفال)

BABY TOOTH EXTRACTION FORCEPS

DG1-1711

English Pattern
Klein Fig. 1
Pedodontic,
Upper Incisors



DG1-1715

English Pattern
Klein Fig. 2
Pedodontic Upper Premolars



DG1-1712

English Pattern
Klein Fig. 3
Pedodontic Upper
Molars Universal



DG1-1716

English Pattern
Klein Fig. 4
Pedodontic Upper root



DG1-1713

English Pattern Klein Fig. 5
Pedodontic Lower Incisors.



DG1-1717

English Pattern
Klein Fig. 6
Pedodontic Lower Molars



DG1-17104

Set of 7 Pieces




DG1-1718

English Pattern
Klein Fig. 7
Pedodontic Lower Molars,



DG1-1719

Comellas Forceps
Anatomical Handle



DG1-1720

Koehler
Luxation forceps
right & left molars



DG1-1721

Fig. 00
Ottavi inferiori
Lower third molars



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)
提取强制 (英语模式)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (نمط اللغة الإنجليزية)

DG1-1721A DEEP GRIPING EXTRACTING FORCEPS SET



• No Fracture • No Tissue Damage • Extraction With Easy



DG1-1722
Fig. 34
Upper Anterior



DG1-1723
Fig. 35
Upper Universal



DG1-1724
Fig. 321
Lower Universal



DG1-1725
Fig. 36A
Lower Premolars
Deep-Gripping



DG1-1726
Fig. 67
Upper Third Molars



DG1-1727
Fig. 79 T
Lower Third Molars

A-TRAUMATIC EXTRACTION FORCEPS
Improved tapered beaks for complete root contact.
Cross serrations for a never slip grasp.
Easily enter the sub gingival.

DG1-1727A Traumatic Extracting Forceps Dual Grip Teeth Serration Set of 6 Pcs



DG1-1728
Fig. 336
lower premolars
deep-gripping



DG1-1729
Fig. 1
Upper Anterior



DG1-1730
Fig. 35
Upper
Universal



DG1-232 D
Fig. 36
Lower Premolars
Deep-Gripping



DG1-1731
Fig. 67
Upper Third
Molars




DG1-1732
Fig. 79A
Apex Forceps Lower
Universal Deep Root Pick

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)
Daviers - forme américaine
Fórceps de extracción - forma americana
提取强制 (美国模式)

Щипці екстракційні (американська форма)
Pinzeper per estrazione - forma americana
ستخراج الملقط (النمط الأمريكي)



ORTHODONTIC INSTRUMENTS

Инструменты ортодонтические
Instrumentos Ortodonticos
Orthodontische Instrumente
正畸仪器

Инструменти ортодонтични
Instruments Orthodontiques
أدوات تقويم الأسنان



DG1-1752
Height Gauge 0.18
Aluminum



DG1-1753
Height Gauge 0.22
Aluminum



DG1-1754
Moveable
Height Gauge



DG1-1755
Bone Bracket Gauge



DG1-1756
Endo Finger Ruler
Span Measure Scale



DG1-1757
Bone Thickness Gauge
40mm, 18cm



DG1-1758
Osteometer
Measuring Forceps For Jaw
Ridge Width measurement, scale
on both sides, rustfree, autoclavable



DG1-1759
Fig. 1
Iwanson
spring caliper
For metal and
porcelain.



DG1-1760
Fig. 2
Iwanson
spring caliper
For wax.



DG1-1761
Castroviejo Calipers
20mm, 45°
9cm



DG1-1762
Castroviejo Calipers
20mm
9cm



A



B



C



D

DG1-1763 ARCH
Generic 1Pc Dental Orthodontic Arch
Forming Turret 0 Degree Without Torque
available color



DG1-1764
Direct Bond Bracket
Holder Tweezers



DG1-1765
Buccal tube tweezer
(Extra Fine Tip)



DG1-1766
Direct Bonding
Bracket Holders
13.5cm



DG1-1767
Direct Bonding
Bracket Holders
13.6cm



DG1-1768
Direct Bonding
Bracket Holders
13.5cm



DG1-1769
Deluxe Bracket
Holder



Инструменты ортодонтические
Instrumentos Ortodonticos
Orthodontische Instrumente
正畸仪器

Інструменти ортодонтичні
Instruments Orthodontiques
أدوات تقويم الأسنان

Available in Hollow Handle Ø9.5mm Handle



DG1-1770 Ligature Director/Scaler



DG1-1771 Dull Scaler/Band Pusher



DG1-1772 Ligature Director/Band Pusher



DG1-1773 Ligature Director Ld1



DG1-1774 Ligature Director/Plugger LDP




DG1-1775 Ligature Director Tucker LDT



DG1-1776 Twiel-on elastic module instrument



DG1-1777 .030" and .018" Distal Bender



DG1-1782 SET

Wamkey Set of 3 Pcs Including Drill



DG1-1778
Band seater
auto cleavable
Triangle Shape



DG1-1779
Band seater
auto cleavable
Square Shape



DG1-1780
Band seater
auto cleavable
Round Shape




DG1-1781
Mershon Band
Pusher
14cm

Щипцы для ортодонтии и протезирования
Pincers pour l'orthodontie et la prothèse
Zangen für Kieferorthopädie und Prothetik
牙齒矯正鉗


Щипці для ортодонції та протезування
Alicates para ortodoncia y prótesis
كماشة لتقويم الأسنان بروتستيس



DG1-1783
Aderer
Large wire bending
and loop forming pliers
13cm



DG1-1784
Angle Plier
Fig. 1
13.5cm



DG1-1785
Adams Plier
Adams pliers for bending wire
arches and bends right angles
for wires up to .022 inch, .025 tip
to beak. Length 13cm



DG1-1786
How Plier
Universally usable tool for the
multiband technology for wires
up to .028 inch. Length 14.5cm



DG1-1787
How Plier curved
Universally usable tool for
the multiband technology for
wires up to .028 inch. Length 14.5cm



DG1-1788
Weingart Plier medium
Weingart plier that are ideal due to
their curved tip for setting sheets and
is capable of bending the ends.
Length 14.5cm



DG1-1789
Tweed Loop
forming Pliers
14cm



DG1-1790
Nance Loop
forming Pliers
14cm



DG1-1791
Band Removing Plier
The curved blade of the band
removing plier can strip easily
and gently remove directly to
the tooth. Length 14cm



DG1-1792
Bracket Removing Plier curved for
easy removal of brackets in the
posterior region. Length 14.5cm



DG1-1793
Bracket removing plier for
comfortable and quick
removal of the brackets.
Length 13.5cm



DG1-1794
Mini Ligature cutters T.C
fine tip 13cm



DG1-1795
Large
Fluch-cut distal and cutting
Pliers T.C with safety hold
fitted with un-breakable,
corrosion resistant and
durable stainless steel tips 13cm



DG1-1796
Small
Fluch-cut distal and cutting
Pliers T.C with safety hold
fitted with un-breakable,
corrosion resistant and
durable stainless steel tips 13cm

Щипцы для ортодонтии и протезирования
Pincettes pour l'orthodontie et la prothèse
Zangen für Kieferorthopädie und Prothetik
牙齒矯正鉗

Щипці для ортодонції та протезування
Alicates para ortodoncia y prótesis
كماشة لتقويم الأسنان بروتستيس



DG1-1797
Jarabak Pliers



DG1-1798
Light Wire Pliers With Groove



DG1-1799
Lingual Arch Forming Pliers



DG1-1800
Hammerhead
NiTi Distal Cinch Back Pliers



DG1-1801
Bird Beak Pliers



DG1-1802
Bracket Removing Plier
Angled Tc - Inserts



DG1-1803
Dental Orthodontic Pliers
Crimpable Hook Placement



DG1-1804
Crimpable Hook Placement



DG1-1805
Rectangular Arch
Forming Pliers




DG1-1806
Direct angled band plier
anterior posterior bracket
removing plier



DG1-1807
Wald sachs Wire Bending
and Cutting Pliers
15cm



DG1-1808
Universal Dolphin Style
Wire Bending & cutting pliers



DG1-1809
Arrow Clasp Forming Pliers
12.5cm



DG1-1810
Schwarz
13cm



DG1-1811
Band Removing Plier
The curved blade of the band
removing plier can strip easily
and gently remove directly to
the tooth. Length 14cm



DG1-1812
Orthodontic Set of 21 Pcs



DG1-1813
Orthodontic Set of 18 Pcs



DG1-1814
Orthodontic Set of 10 Pcs

Щипцы для ортодонтии и протезирования
Pincettes pour l'orthodontie et la prothèse
Zangen für Kieferorthopädie und Prothetik
牙齒矯正鉗

Щипці для ортодонції та протезування
Alicates para ortodoncia y prótesis
کماشة لتقويم الأسنان بروتستیس

THERMO PLIERS & MATHIEU NEEDLE HOLDERS



DG1-1815
Hole Punch Plier
Creates half-moon cutouts in clear aligners to allow space for bonded buttons used for elastics



DG1-1816
Vertical Punch Plier
This versatile plier emphasizes on over-correction of rotations by creating an indentation at the mesial or distal line angles on the facial and/or lingual of a specific tooth. Can also be used to create a rotational couple using vertical indents on facial and lingual surfaces on opposite sides of a tooth.



DG1-1817
Horizontal Punch Plier
Designed to produce an indentation to accent individual root torque. Can also be used to enhance retention of clear aligner or retainers by placing an indent into a tooth's undercut or below a bonded attachment.



DG1-1818
Tear Drop Punch Plier
Designed to create a notch at the gingival margin of clear aligners anywhere elastic hooks are required within the arch, making sure the elastic stays in place once the aligner is seated.



DG1-1819
Dimple Punch Plier
Creates a tiny rounded dimple for spot contacts to accent individual tooth movements or added aligner retention. Typically used to accent movement on any tooth that requires a slight "nudge"



DG1-1820
Eraser Plier
The eraser plier is designed to remove and flatten any unwanted indentations placed on clear aligners by other punch pliers.



DG1-1821
Groove thermal forming plier
Designed for making grooves on clear appliances, to hold the elastic, prevent it falling off



DG1-1822 Yellow
Thermoplier Micro-Ramp Pliers
Tip, torque or rotate smaller teeth laterally; Allow contact close to incisal edge. On the correct location on the target tooth



DG1-1823 Pink
Thermoplier Small Circle Pliers
Increase the size (depth) of small existing force point. Location of the existing force point



DG1-1824
Crown Beebe Scissors
4.5 Str.



DG1-1825
Mathieu Standard Tip



DG1-1826
Elastomeric ligature placing forceps with groove 14cm



DG1-1827
Mathieu Narrow Tip



DG1-1828
Mathieu Smaha Ultra fine tip



DG1-1831
Str.




DG1-1832
Cvd.



DG1-1829
Mathieu Standard Tip Tungsten Carbide



DG1-1830
Mathieu Narrow Tip Tungsten Carbide




Smaha Elastomeric ligature placing forceps 140 mm

Инструменты для имплантологии
Implantologie
Implantologie
种植牙

Інструменти для імплантології
Implantoloia
زراعة الأسنان

DG1-1832A SINUS LIFT ELEVATOR & MEMBRANE LIFT TIPS
Available in Hollow Handle Ø9.5mm



With curette the sinus membrane is further detached over a broad area in the inner region, especially along the roof of the palate.

With curette the sinus membrane is further elevated in the upper and posterior part of the maxillary sinus wall.

The "Schneider's Membrane" of the sinus can be detached with curette by 1–2mm along the total cortical bone.

Curette is particularly useful in the lower part for detaching the membrane further. The same goes for the caudal area, especially when the membrane is 4–5 mm above the base of the maxillary sinus and the maxillary sinus is built up with granular material.

Curette makes it possible to lift the membrane further in mesial region. One working tip is especially designed for the right maxillary sinus and the other for the left maxillary sinus.

Инструменты для имплантологии
Implantologie
Implantologie
种植牙

Інструменти для імплантології
Implantoloia
زراعة الأسنان

AAY8-522 SINUS LIFT INSTRUMENTS SET OF 10 PCS
Available in Hollow Handle Ø9.5mm



With curette the sinus membrane is further detached over a broad area in the inner region, especially along the roof of the palate.

With curette the sinus membrane is further elevated in the upper and posterior part of the maxillary sinus wall.

The "Schneider's Membrane" of the sinus can be detached with curette by 1–2mm along the total cortical bone.

Curette is particularly useful in the lower part for detaching the membrane further. The same goes for the caudal area, especially when the membrane is 4–5 mm above the base of the maxillary sinus and the maxillary sinus is built up with granular material.

Curette makes it possible to lift the membrane further in mesial region. One working tip is especially designed for the right maxillary sinus and the other for the left maxillary sinus.

Инструменты для имплантологии
Implantologie
Implantologie
种植牙

Инструменти для імплантології
Implantologia
زراعة الأسنان

DG1-1855 Vista Tunneling Set of 8 Pcs
Available in 8.0mm - 9.5mm Hollow Handle



DG1-1856 VISTA1



DG1-1857 VISTA2




DG1-1858 VISTA3



DG1-1859 VISTA4



DG1-1860 VISTA5



DG1-1861 VISTA6



DG1-1862 VISTA7

TUNNELING GUM LIFT EXTRA THIN, NOT FOR EXCESSIVE PRESSURE, Micro Serrated




DG1-1863 PP Buser 3.8 - 4.3mm



DG1-1864 Atraumatic Extraction Periotomes Set



DG1-1865 Periotome Pt1



DG1-1866 Periotome Pt2




DG1-1867 PPAEL 1.7mm Curved



DG1-1868 PPAELA 2mm



DG1-1869 PT-4 Single Ended Straight



DG1-1870 PT-5 Single Ended Curved

PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie
牙周仪器

Інструменти пародонтологічні
instrumentos para periodoncia
الصكوك اللثة

Available in Hollow Handle Ø9.5mm

DG1-1871 SET **PERIODONTAL CHISEL SET OF 6 PCS**



DG1-1872 TG



DG1-1873 13K/TG Titanium Plasma Coated Point



DG1-1874 13K/13KL



DG1-1875 FEDI3



DG1-1876 KK 15/16 Kirkland Knife



DG1-1877 Orban Knife 1/2

DG1-1878 SET **PERIOSTEAL ELEVATORS SET OF 7 PCS**



DG1-1879 Molt



DG1-1880 Molt



DG1-1881 Molt



DG1-1882 Molt



DG1-1883 P 24 G



DG1-1884 Molt N. 9



DG1-1885 Periosteal Pichard Pr3

Кюрети
Curetas
Küretten
水果

Кюрети
Curette
كرتس

Available in Hollow Handle Ø8.0mm - 9.5mm

DG1-1886 SET LUCAS CURETTES SET OF 5 PCS



DG1-1892 SET LUCAS CURETTES SERRATED SET OF 5 PCS



Инструменты для имплантологии
Implantologie
Implantologie
种植学

Інструменти для імплантології
Implantologia
زراعة الأسنان

DG1-1898 BONE GRAFTING SET

Blue Titanium Coated TIPS

Available in Hollow Handle Ø9.5mm



DG1-1899

Bone Graft 1.3 / 1.7 mm Titanium Coated tips 8. 11. 13. 15. 18 mm

DG1-1900

Bone Graft 1.5 / 2.0 mm Titanium Coated tips 8. 11. 13. 15. 18 mm

DG1-1901

Bone Graft 2.4 / 2.7 mm Titanium Coated tips 8. 11. 13. 15. 18 mm

DG1-1902

Bone Graft 3.3 / 4.0 mm Titanium Coated tips 8. 11. 13. 15. 18 mm

DG1-1903

Bone Graft Titanium Coated tips

DG1-1904

Bone Graft / Carrier Titanium Coated tips

DG1-1905 Gingival Retraction Cord Packers Serrated Set of 5 Pcs

Available in Hollow Handle Ø9.5mm

DG1-1906

Gingival Cord Packer Serrated

DG1-1907

Gingival Cord Packer Serrated GCPCSI 1

DG1-1908

Gingival Cord Packer Serrated GCP 113

DG1-1909

Gingival Cord Packer 2.0mm GCP 171

DG1-1910

Retraction Cord Packer 3.5mm GCP 170

Инструменты для имплантологии
Implantologie
Implantologie
种植牙

Інструменти для імплантології
Implantologia
زراعة الأسنان

DG1-1911 SET Tunnelling Knife Set of 4 Pcs



DG1-1912 Tunnelling Knife 1



DG1-1913 Tunnelling Knife 2



DG1-1914 Tunnelling Knife 3



DG1-1915 Tunnelling Knife 4

DG1-1916 Sinus Instruments Set of 4 Pcs



XSE1L



XSE2L




XSE3L



XSE4L



DG1-1917
Soft Brushing Kit
Set of 6 Pcs



DG1-1918 SET Flexible Perioste Set of 6 Pcs

Инструменты для имплантологии
Implantologie
Implantologia
种植牙

Інструменти для імплантології
Implantologia
زراعة الأسنان

CONVEXES



2,7 Ømm, 3,2 Ømm, 3,7 Ømm, 4,2 Ømm, 5 Ømm

DG1-1919
SET OF 5 PIECES
Convexes

CONCAVES



2,7 Ømm, 3,2 Ømm, 3,7 Ømm, 4,2 Ømm, 5 Ømm

DG1-1920
SET OF 5 PIECES
Concaves

CONIQUES



2,7 Ømm, 3,2 Ømm, 3,7 Ømm, 4,2 Ømm, 5 Ømm

DG1-1921
SET OF 5 PIECES
Coniques



DG1-1922 SET
Osteotome Convexes
Straight Curved
SET OF 11 PIECES



DG1-1923 SET
Osteotome Concaves
Straight Curved
SET OF 11 PIECES



DG1-1924 SET
Osteotome Coniques
Straight Curved
SET OF 11 PIECES

CONVEXES



2,7 Ømm, 3,2 Ømm, 3,7 Ømm, 4,2 Ømm, 5 Ømm

DG1-1925
SET OF 5 PIECES
Convexes Curved

CONCAVES



2,7 Ømm, 3,2 Ømm, 3,7 Ømm, 4,2 Ømm, 5 Ømm

DG1-1926
SET OF 5 PIECES
Concaves Curved

CONIQUES



2,7 Ømm, 3,2 Ømm, 3,7 Ømm, 4,2 Ømm, 5 Ømm

DG1-1927
SET OF 5 PIECES
Coniques Curved

Инструменты для имплантологии
Implantologie
Implantologie
种植牙

Инструменти для імплантології
Implantologia
زراعة الأسنان

BONE MILL BONE WELL



DG1-1928
Bone Crusher




DG1-1929
Bone Mill



DG1-1930
Basins For Bone



DG1-1931
Gold Mill (Small)



DG1-1932
Gold mill (Large)

TORQUE WRENCH & RATCHET



DG1-1933 Star Torque Wrench 6.35mm Universal



DG1-1934 Round Torque Wrench 4.0mm Universal



DG1-1935 Hex Driver With 4.0mm Driver




DG1-1936 Star Torque Wrench 6.35mm Universal


TISSUE PUNCH STRAIGHT & CURVED



Tissue Punch
Straight



Tissue Punch
Curved



BONE SCRAPER STRAIGHT & CURVED



DG1-1949
Bone Scraper
Straight



DG1-1950
Bone Scraper
Curved

BONE INJECTORS, COLLECTORS & SUCTION TUBE



DG1-1951 Star Ratchet Wrench 6.35mm Universal



DG1-1952 Star Ratchet Wrench 6.35mm Universal



DG1-1953 Star Ratchet Wrench 4.0mm Universal



DG1-1954 Round Ratchet Wrench 3.35mm Universal



DG1-1955
Implant Depth Gauge



DG1-1956
Implant Depth Gauge




DG1-1957 Double Ended Osteotomy Depth Gauge



DG1-1958 Double Ended Osteotomy Depth Gauge



DG1-1959
Ø 2.5mm
16cm
I.D 2.5mm
O.D 3.5



DG1-1960
3.5mm
16cm
I.D 4.5mm
O.D 3.5


Инструменты для имплантологии
Implantologie
Implantologia
زراعة الأسنان
种植牙学

Инструменти для імплантології
Implantologia
زراعة الأسنان



DG1-1961 Suction Tubes 18cm Blue Titanium colour
DG1-1962 Suction Tubes 18cm Blue Titanium colour
DG1-1963 Suction Tubes 18cm Blue Titanium colour
DG1-1964 Suction Tubes 18cm Blue Titanium colour
DG1-1965 Suction Tubes 18cm Blue Titanium colour
DG1-1966 Suction Tubes 18cm
DG1-1967 Suction Tubes 18cm
DG1-1968 Suction Tubes 18cm
DG1-1969 Suction Tubes 18cm
DG1-1970 Suction Tubes 18cm

Frazier Suction Tube Blue Titanium



DG1-1971 Ø 3 mm length 18cm Titanium color
DG1-1972 Ø 4 mm length 18cm Titanium color
DG1-1973 Ø 5 mm length 18cm Titanium color

DG1-1974 Ø 3 mm length 18cm
DG1-1975 Ø 4 mm length 18cm
DG1-1976 Ø 5 mm length 18cm

DG1-1977 Fig. 6FR Frazier 15cm - Ø 2 mm

DG1-1978 fine Ø 9mm Suction tubes to collect the extracted bone material 180mm

DG1-1979 17.5cm Ø 20mm

DG1-1980 18cm Ø 27mm

DG1-1981 coarse Ø 5mm

Инструменты для имплантологии
Implantologie
Implantologia
زراعة الأسنان
种植牙

Инструменти для імплантології
Implantologia
زراعة الأسنان

PERIODONTAL SURGERY KIT



Available in Hollow Handle 8mm, Ø9.5mm



Periodontal Chisels 4mm / 5mm



Curette curved,
Langer 1/2



Curette straight,
Langer 5/6



Hygienist H6/H7 Scaler



Periosteal Prichard PR3




Rhodes Back Action
C36/37



Ochsenbein Bone Chisel



Lucas Bone Curette,
84



Periodontal Probe,
UNC 15
CNC North Carolina,
graduation 1-15
in mm steps



Round Scalpel Handle [8mm Diameter, 14.5cm Length]
Scalpel Handle straight, for blades # 10 -15, 14.5 cm



Micro Surgical Scalpel Handle, titanium, 13.5 cm



Iris Scissors T.c 11.5cm



Cooley
18 cm




Noyes Micro Scissors CVD 17cm



Castroviejo TC 18cm



DG1-1982
PERIODONTAL SURGERY KIT
of 17 Pcs



Инструменты для имплантологии
Implantologie
Implantologia
زراعة الأسنان
种植学

Інструменти для імплантології
Implantologia
زراعة الأسنان

DENTAL IMPLANT SURGERY KIT

DENTAL IMPLANT GUIDED BONE REGENERATION GBR KIT



Graft apacker, short lexible, Strai, titanium, Serrated 3mm/4mm



Graft Packer, Long, Titanium 8. 11. 13. 15. 18 mm Smooth 3.3 / 4.0 mm



Buser Periosteal Elevator, Flat Tip, Titanium



Kirkland Periodontal Chisel, Titanium 13KL/LR



Woodson Periosteal Elevator, Titanium



Prichard Pr3 Periosteal Elevator, Titanium



Bone Graft 3.3mm/ 4.0mm Titanium Coated Tips



Scalpel Handle straight, #3



Micro Surgical Scalpel Handle



Fig. 3 12.5cm



Goldman-Fox 13 cm Saw edge, Curved



Forceps 18cm



Bone Scraper Straight



Castroviejo Needle Holder TC Straight 18cm



DG1-1983 Set of 15 Pcs with Cassette Tray



Basins For Bone


Инструменты для имплантологии
Implantologie
Implantologie
种植牙

Інструменти для імплантології
Implantologia
زراعة الأسنان


DENTAL IMPLANT SURGERY KIT DENTAL IMPLANT SURGERY KIT



CP-12 University of North Carolina



Columbia 4R/4L Posterior Application




Orban knives 1/2



24G 4mm 4mm



Molt N. 9 4mm 7mm



Lucas Curette 3.5mm



Rhodes back-Action 36/37



3.5mm

13K/13KL

4.5mm



PR-3 Prichard



DG1-1984
Set of 19 Pcs
with Cassette Tray



Iris Scissors T.c
11.5cm



Mayo-Hegar
16cm



Castroviejo Needle
Holder TC
Straight 18cm



Castroviejo
Scissors
TC Curved 18cm



Forceps
18cm



Minnesota
Large
16cm



Scalpel Handle
straight, for blades



Locks At 90°, 60°, 30°,
Straight, 120°, 150°, 180°
Angles



Micro Surgical
Scalpel Handle



Scalpel Handle
Adjustable

Инструменты для имплантологии
Implantologie
Implantologie
种植牙学

Інструменти для імплантології
Implantologia
زراعة الأسنان

DENTAL IMPLANT SURGERY KIT




Castroviejo Needle
Holder TC
Straight 18cm

Castroviejo
Scissors
TC Curved 18cm



DG1-1985 set of 11 Pcs
With Cassette, 114×202×36mm20



Tunnelling Knife 1

Tunnelling Knife 2

Tunnelling Knife 3

Tunnelling Knife 4



Micro Surgical
Scalpel Handle




Scalpel Handle
straight, #3



Corn Suture Plier,
cvd. 16 cm



Crile-Wood Needle
Holder TC 15cm



DG1-1986 PRF/GRF INSTRUMENTS SET

Bone Grafting Applicator, 12mm/16mm, DE

Bone Grafting Applicator, 6mm/10mm, DE

Bone Compactor Ø6.0mm/8.0mm, DE

Bone Compactor Ø3.0mm/4.0mm, DE



PRF Box, with lid and Cover



Tube Holder, for 12 tubes



Mixing Bowl for bone material,
28mm, Ø65mm



Goldman-Fox, Str. TC Scissor, 12.5cm



Tweezer Forceps T.C



Splinter T.C Angular Tweezer Forceps



DG1-373
Bone Chip Maker 2

Mini Tissue Punch & Trephine Bur & Saw Set



DG1-374
Tissue Punch & Trephine Bur & Saw Set



DG1-375
Dental Implant Guide System



DG1-376
Bone & Ridge Expansion Kit

Инструменты для имплантологии
Implantologie
Implantologie

Инструменти для імплантології
Implantologia
زراعة الأسنان



DG1-381

Trimmer Kit

The Kit Includes:

1) Trephine Drills

3.0mm, 4.0mm
5.0mm, 6.0mm
7.0mm

2) Tissue Punch

3.5mm, 4.0mm
4.5mm, 5.0mm
5.5mm

3) Guide Punch

4.0mm, 5.0mm

4) Saw Disks

7.0mm, 10.0mm

5) Sterilization Cassette.



DG1-382

Bone Compression Kit

The Kit Includes:

Convex Compression Screws

1.6/2.7mm, 2.0/3.3mm
2.4/3.6mm, 3.1/4.1mm

Cancave Compression Screws

3.0mm, 3.5mm, 4.0mm

Spade Marking Drill

Motor Mount Drivers Long & Short
Hand Driver
Ratchet Adapter
Ratchet Wrench



DG1-383

Surgical Drill Kit

The Kit Includes:

1) Twist Drills

2.0mm, 2.8mm
3.2mm, 3.65mm
4.2mm, 4.8mm, 5.2mm

With Laser Depth Markings on 6.0mm,
8.0mm, 10.0mm, 11.5mm, 13.0mm & 16.0mm

2) Ratchet / Wrench Hex Drivers (Color Coded)

1.25mm Short / Long, 2.42mm Short / Long

3) Parallel Pins

Long & Short

4) Motor Mount Hex

1.25mm, 2.5mm

5) Ratchet Wrench 4mm

6) Sterilization Cassette / Box / Organizer

Инструменты для имплантологии
Implantologie
Implantologie

Инструменти для імплантології
Implantologia
زراعة الأسنان



DG1-384

Universal Kit 50 PCS

The Kit Includes:

16pcs Conical Drills with Stoppers

1.8mm / 2.2mm - S1, 1.8mm / 2.2mm - S3
1.8mm / 2.2mm - S5, 1.8mm / 2.2mm - S7
2.2mm / 2.7mm - S1, 2.2mm / 2.7mm - S3
2.2mm / 2.7mm - S5, 2.2mm / 2.7mm - S7
2.3mm / 3.3mm - S1, 2.3mm / 3.3mm - S3
2.3mm / 3.3mm - S5, 2.3mm / 3.3mm - S7
2.8mm / 4.5mm - S1, 2.8mm / 4.5mm - S3
2.8mm / 4.5mm - S5, 2.8mm / 4.5mm - S7
S = Stopper Length

Bone Expander Screws

1.6mm / 2.6mm, 1.9mm / 3.0mm
2.9mm / 3.4mm, 2.3mm / 3.8mm
2.9mm / 4.3mm,

Titanium Coated Drills

Drill Bit 2.0mm L 18.50, Drill Bit 2.8mm L 18.50
Drill Bit 3.65mm L 18.50, Drill Bit 4.2mm L 18.50
Drill Bit 4.5mm L 18.50, Drill Bit 4.8mm L 18.50
Drill Bit 5.2mm L 18.50

8pcs Trephine Drills

Int: 2.0mm, Ext: 3.0mm, Int: 3.0mm, Ext: 4.0mm
Int: 4.0mm, Ext: 5.0mm, Int: 5.0mm, Ext: 6.0mm
Int: 6.0mm, Ext: 7.0mm, Int: 7.0mm, Ext: 8.0mm
Int: 8.0mm, Ext: 9.0mm, Int: 9.0mm, Ext: 10.0mm
Compression Screw 2.8 / 4.0mm, Compression Screw 4.0mm
Driver 2.5mm
Compression Screw 2.0 / 2.8mm, Compression Screw 2.8 / 3.4mm
Tissue Punch 4.5mm, Tissue Punch 5.8mm
Pilot Drill 2.0mm
Drill Extender
Ratchet Driver 2.5mm, Hand Driver 2.5mm
Torque Wrench 6.35mm (10-40Ncm)
Depth Gauge
Sterilization Cassette / Box



DG1-385

Conical Implant Drill Kit

The Kit Includes:

Drills Sizes:

1.8mm / 2.2mm, 2.2mm / 2.7mm
2.3mm / 3.3mm, 2.7mm / 3.7mm
2.7mm / 4.5mm

Stopper Sizes:

1.0mm / L 13.0, 2.0mm / L 11.5
3.0mm / L 10.0, 5.0mm / L 8.0
7.0mm / L 6.0

Overall Drills Length:

13.0mm (32mm with Shaft)

Инструменты для имплантологии
Implantologie
Implantologie

Инструменти для імплантології
Implantologia
زراعة الأسنان



DG1-377

Compression Expander Screws Kit

The Kit Includes:

Kit Includes

3x Bone Compression Screws

1x Sinus Elevation Screw (3.5mm)

Screws are compatible with any brands 4.0mm Torque Wrench / Ratchet Wrench / Adapter's

Kit Advantages

Screws are Laser marked to indicate the depth of 8mm, 10mm & 13mm



DG1-378

Hex Drivers Adapter Kit

The Kit Includes:

4 Pcs Hex Drivers Kit Fits Ratchet / Torque Wrench

2 x 25mm & 2x 2.5mm Long & Short Driver

Manufactured from High Grade Stainless Steel - Highly Polished



DG1-379

Bone Expander Kit

The Kit Includes:

1) Expander Screws

1.6 / 2.6 mm, 1.9 / 3.0 mm, 2.1 / 3.4 mm

2.6 / 3.8 mm, 2.9 / 4.3 mm

With Laser Depth Markings on 8mm, 10mm & 11.5mm.

2) Saw Disks

7.0 mm, 10.0 mm

3) Pilot Drill 1.6mm


4) Hand piece Adapter

5) Ratchet Extender

6) Hand Driver

7) Ratchet Wrench 4.0mm

8) Sterilization Cassette / Box.



DG1-380

Bone Compact Kit

The Kit Includes:

1) Twist Drills

2.0mm, 2.5mm, 2.8mm 3.2mm,

3.65mm, 4.2mm, 5.2mm

With Laser Depth Markings

on 6.0mm, 8.0mm, 10.0mm,

11.5mm, 13.0mm & 16.0mm

2) Ratchet / Wrench Hex Drivers

(Color Coded)

1.25mm Short / Long, 2.42mm Short / Long

3) Contra Angle Hex Drivers

1.25mm, 2.42mm

4) Parallel Pins Long & Short

5) Drill Extender

6) Marking Drill

7) Ratchet Wrench 6.35mm

8) Sterilization Cassette / Box / Organizer

Инструменты для имплантологии
Implantologie
Implantologia

Інструменти для імплантології
Implantologia
زراعة الأسنان

BONE EXPANDER SCREW & SURGICAL KITS



DG1-386

Bone Universal Kit
Drill: 2mm, 2.8mm, 3.2mm, 3.65mm, 4.2mm, 4.5mm, 4.8mm, 5.2mm
Trephines Bit: 3.0mm, 4.0mm, 5.0mm, 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm.
Driver: 2.5mm, 2.5mm Hand Pcs
1.25mm, 1.25mm Hand Pcs
Short Driver 2.5mm
Long Driver 2.5mm
1.25mm, Short Hex Driver Hand Pcs
1.25mm, Long Hex Driver Hand Pcs
Torque Wrench 4.0mm



DG1-387

Bone Expander Kit
Expander: 2.6mm, 3.0mm, 3.4mm, 3.8mm, 4.3mm
Ratchet Driver
Parallel Pin
Pilot Drill
Drill 1.4mm
Drill Driver 4mm
Hand Driver 4mm
Ratchet 4.0mm



DG1-388

Bone Universal Kit
Drill: 2mm, 2.8mm, 3.2mm, 3.65mm, 4.2mm, 4.5mm, 4.8mm, 5.2mm
Trephines Bit: 3.0mm, 4.0mm, 5.0mm, 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm.

Bone Compression: 2mm, 3.3mm, 2.8mm, 4.2mm, 3mm, 3.5mm
Driver: 2.5mm, 2.5mm Hand Pcs
1.25mm, 1.25mm Hand Pcs
Short Driver 2.5mm
Long Driver 2.5mm
1.25mm, Short Hex Driver Hand Pcs
1.25mm, Long Hex Driver Hand Pcs
Torque Wrench 6.35mm 1 Pcs



DG1-389

Conical Drills: 20 Pcs
Tissue Punches: 3.0mm, 3.50mm, 4.0mm, 4.5mm, 5.0mm, 6.0mm, 6.50mm
Bone Expander: 1.6/3.8mm, 1.9/3.0mm, 2.1/3.6mm, 2.6/3.8mm, 2.9/4.3mm
Drill Adapter: 1 Pcs
1.25mm, 2.50mm Hex Driver Hand Pcs

Pilot Drill: 1.6mm
Drill Bit: 1.4mm, 1.6mm, 1.8mm, 2.0mm, 2.2mm
Ratchet Driver: 4mm/6.35, 1.25/2.50
Rimer: 3.2/4.7mm, 3.7/4.7mm
Parallel Pin: 4.5mm, 2.5mm
Double step Drill: 2.4/2.8mm
Depth Gauge: 1 Pcs
Torque Wrench 6.35mm Universal: 1 Pcs

Инструменты для имплантологии
Implantologie
Implantologie

Инструменти для імплантології
Implantologia
زراعة الأسنان



DG1-394
Drills kit (Black Coated)
The Kit Includes:
Dental Drills Implant Kit 8 Pcs
Drill Sizes:
2.0mm, 2.5mm, 2.8mm, 3.2mm,
3.65mm, 4.2mm, 5.2mm, 5.5mm
Marking Depths on the drills
for accurate drilling,
Sizes 6.0mm, 8.0mm, 10.0mm,
11.5mm, 13.0mm & 16.0mm.



DG1-395
Drill Kit (Polished)
The Kit Includes:
Dental Drills Implant Kit 8 Pcs
Drill Sizes:
2.0mm, 2.5mm, 2.8mm, 3.2mm,
3.65mm, 4.2mm, 5.2mm, 5.5mm
Marking Depths on the drills
for accurate drilling,
Sizes 6.0mm, 8.0mm, 10.0mm,
11.5mm, 13.0mm & 16.0mm.

The Kit Includes:

- 1.8mm Pilot Drill
- 1525 (2.0mm)
- 2535 (3.0mm)
- 3545 (4.0mm)
- 4555 (5.0mm)
- 1828 (2.3mm)
- 2838 (3.3mm)
- 3848 (4.3mm)
- 4858 (5.3mm)
- 2228 (2.5mm)
- 3238 (3.5mm)
- 4248 (4.5mm)
- 5258 (5.5mm)



DG1-396
Osseodensification Burs Kits

Инструменты для имплантологии
Implantologie
Implantologie

Инструменти для імплантології
Implantologia
زراعة الأسنان



DG1-397
Osseodensification Burs Kit
The Kit Includes:
12pcs Conical Drills



DG1-398
Osseodensification Burs Kit
The Kit Includes:
18pcs Conical Drills Kit



DG1-399
Osseodensification Burs Kit
The Kit Includes:
13pcs Conical Drills



DG1-400
Positioning Kit
The Kit Includes:
11pcs Drills Guide & Implants



DG1-401
Premium Surgical Drill Kit
The Kit Includes:
44pcs Drills



DG1-402
Premium Surgical Drill Kit
The Kit Includes:
51pcs Drills



DG1-403
Premium Surgical Drill Kit
The Kit Includes:
26pcs Drills



DG1-404
Premium Surgical Drill Kit
The Kit Includes:
24pcs Drills



DG1-405
Bone Expander and Chisel Kit
The Kit Includes:
16pcs Chisel and Expander



DG1-406
Sinus and Stopper Kit
The Kit Includes:
14pcs



DG1-407
Sinus and Stopper Kit
The Kit Includes:
14pcs



DG1-408
Sinus and Stopper Kit
The Kit Includes:
14pcs



DG1-408A
Sinus and Stopper Kit
The Kit Includes:
16pcs



DG1-409
Dental Implant Titanium Bushes Kit
The Kit Includes:
6pcs Bushes



DG1-410
Lance Drill Kit
The Kit Includes:
4pcs Drill



DG1-411
Bone Screw Kit
The Kit Includes:
5pcs Screw



DG1-412
Screwdrivers Kit
The Kit Includes:
8pcs Screwdriver



DG1-413
Saw And Trephine Kit
The Kit Includes:
6pcs Trephine



DG1-414
Saw And Trephine Kit
The Kit Includes:
6pcs Trephine




DG1-415
Fixture Remover Screw Kit
The Kit Includes:
3pcs Screw



DG1-416
Fixture Remover Screw Kit
The Kit Includes:
3pcs Screw



DG1-417
Lance Drills Kit
The Kit Includes:
3pcs Drills



DG1-418
Bone Chip Maker II With Stopper
The Kit Includes:
10pcs



DG1-419
Dental Residual Root Drill Kit
The Kit Includes:
3pcs Drills



DG1-420
Para Drill Aid
The Kit Includes:
8pcs Drills



DG1-421
Dental Implant Drills



DG1-422
Dental Implant Drills



DG1-423
Dental Implant Drills



DG1-424
Dental Implant Drills



DG1-425
Dental Implant Drills



DG1-426
Dental Guide Positioning Kit



DG1-427
Trephine Drill Kits



DG1-428
Trephine Drill Kits



DG1-429
Trephine Drill Kits



DG1-430
Trephine Drill Kits



DG1-431
Tissue Punch Kit



DG1-432
Tissue Punch Kit



DG1-433
Tissue Punch Kit



DG1-434
Tissue Punch Kit



DG1-435
Bone Saw Kit



DG1-436
Lateral Approach Sinus Lift Kit



DG1-437
Lateral Approach Sinus Lift Drills



DG1-438
Surgical Drill Kit



DG1-439
Universal Densah Bur Kit




DG1-440
Universal Prosthetic Kit



DG1-441
Bone Expander Kit



DG1-442
lateral sinus lift kit





DG1-443
Ruler Measuring gauge



DG1-444
Ruler Measuring gauge



DG1-445
Measuring gauge




DG1-446
Measuring gauge



DG1-447
Measuring gauge



DG1-448
Screws and handle



DG1-449
Screws & Handle Kit





DG1-450
Hammer



DG1-451
Hammer



DG1-452
Hammer



Cement Teeth Filling Plugger Applicator
Carrier Root Canal Endodontic
Material 0.8/1.0/1.2/1.8/2.0/2.8/3.5mm



DG1-453
0.8mm



DG1-454
1mm

IMPRESSION TRAYS


■ Оттисные ложки ■ Відбиткові ложки
 ■ Porte-empreinte ■ Cubetas para impresiones
 ■ Abdrucklöffel ■ صواني الانطباع

UPPER


LOWER



Perforated With Rim Lock DG1-2009 SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L



Perforated With Rim Lock DG1-2010 SET OF 10 PCS
Upper: XS, S, M, L, XL
Lower: XS, S, M, L, XL



DG1-2011 44mm
 DG1-2012 55mm
 DG1-2013 60mm
 DG1-2014 65mm

DG1-2015

DG1-2016 44mm
 DG1-2017 55mm
 DG1-2018 60mm
 DG1-2019 65mm

DG1-2020


DG1-2021

Partial Perforated Set of 5 Pcs. DG1-2022





DG1-2023
182 x 145 x 35 mm
Cassette For 10 Instruments



DG1-2024
182 x 72 x 35 mm
Cassette For 5 Instruments



DG1-2025
Cassette For 7 Instruments



DG1-2026
Cassette Tray 15 Pcs



DG1-2027
Cassette For 20 Instruments



Scalor Tray
DG1-2028 200 x 115 mm
DG1-2029 205 x 150 mm
DG1-2030 290 x 150 mm
DG1-2031 300 x 200 mm



Scalers Trays
DG1-2032 205 x 115 x 17 mm
DG1-2033 205 x 152 x 17 mm
DG1-2034 290 x 115 x 17 mm



Scalers Tray
DG1-2035 200 x 100 mm



DG1-2036
Instruments Tray Standard
Holds up to 7
instruments



DG1-2037
ø mm 80 inox



DG1-2038
mm 200 x 90
Kidney bowls
Bacinella reniforme

DG1-2039
mm 300 x 140
Kidney bowls
Bacinella reniforme



DG1-2040 CL
mm 200 x 90
Cover for kidney bowl
Coperchio per bacinella
reniforme



 DENTOGRADE
MANUFACTURERS

Manufacturers Dental Instruments

Sialkot Pakistan 00923218718711

Email: Info@dentograde.com

Web: <https://dentograde.com>

CE MARK
Approved

