



DENTOGRADE
— MANUFACTURER —

Contact us At:

Phone/WhatsApp: +92 3218718711

+44 7437879455

Email: Dentograde@gmail.com

INSTRUMENT CASSETTES LASER CUT STAINLESS STEEL



- Sizes:- 5 Instruments: 78 mm x 182 mm x 32 mm ...
- 10 Instruments: 142 mm x 182 mm x 32 mm...



COMPOSITE FILLING INSTRUMENTS SILICONE HANDLE / HOLLOW XL HANDLE

SILICONE HANDLE



HOLLOW XL HANDLE




FILLING INSTRUMENTS SILICONE HANDLE



PHOTOGRAPHIC **TITANIUM GLASS** MIRRORS FRONT SURFACE DOUBLE SIDED

- Autoclavable up to 134°
- Front surface mirrors prevent ghosting effect and provides a clear, undistorted view
- Titanium-coated hard surface glass



SCALERS AND CURETTES **SILICONE** HANDLE





PROBES SILICONE HANDLE SINGLE ENDED



47-020-006-SL Probe Silicone Handle Blue Fig 6 ...



47-020-016-SL Probe Silicone Handle Light Green Fig 16 ...



47-020-009-SL Probe Silicone Handle Red Fig 9 ...



47-020-023-SL Probe Silicone Handle Yellow Fig 23 ...

actual size



PERIO PROBES LASER MARKED SILICONE HANDLE SINGLE ENDED



49-181-001-SL Perio Probe Silicone Handle Blue Williams ...



49-181-003-SL Perio Probe Silicone Handle Red 11 ...



49-181-002-SL Perio Probe Silicone Handle Light Green WHO ...



49-181-004-SL Perio Probe Silicone Handle Yellow 12 ...

PERIO PROBES LASER MARKED SILICONE HANDLE DOUBLE END



49-180-001-SL Perio Probe Silicone Handle Blue 23/Williams ...



49-180-002-SL Perio Probe Silicone Handle Light Green 23/WHO ...



49-180-003-SL Perio Probe Silicone Handle Red 23/11 ...



49-180-004-SL Perio Probe Silicone Handle Yellow 23/12 ...

EXCAVATORS SILICONE HANDLE DOUBLE END



49-185-017SL Excavator Silicone Handle Blue No. 17 ...



49-185-019SL Excavator Silicone Handle Light Green No. 19 ...



49-185-018SL Excavator Silicone Handle Red No. 18 ...



49-127-128SL Excavator Silicone Handle Yellow No. 127/128 ...



EXAMINATION INSTRUMENTS

Mirror Handles Stem Silicone Handle
(mirror not included)
47-051-011-SL Blue
47-051-012-SL L.Blue
47-051-013-SL L.Green
47-051-014-SL Yellow
47-051-015-SL Pink

Black Coated Probes / Tweezer
(avoids visual fatigue under light)
49-190-006-BF Probe Fig.6
49-190-009-BF Probe Fig.9
49-190-016-BF Probe Fig.16
49-190-023-BF Probe Fig.23
64-483-150-BF Tweezer College 15cm ...

Mouth Mirrors Front Surface Stem Black Coated
(avoids visual fatigue under light)
47-065-022TB F/S Single Side No.4
47-065-024TB F/S Single Side No.5
47-066-022TB F/S Double Sides No.4
47-066-024TB F/S Double Sides No.5
HH-80-006-BF Hollow Mirror Handle
(mirror not included)



HAHNENKRATT RELAX FRONT SURFACE PLASTIC MIRRORS AUTOCLAVABLE

RELAX - designed with a perfect blend of technology and ergonomics.

- Comfortably light weight for relaxed dental operations
- Comfortable to hold due to the ergonomic design with grip recess
- Keeps quality colour stable, resistant against acids as well as against plaque-indicators
- Comfortable for patients due to the smooth finish
- Hygienic due to the plane transition from mirror surface to mirror casing without any gap
- Silk-mat-finished reflection free surface

Relax Mirrors, 10 pcs/pack
No.4 (22 mm)
7302 Pastel Yellow
7303 Pastel Blue
7304 Pastel Mint
7305 Pastel Lavender ...
7307 Pastel Pink

Relax Mirrors, 10 pcs/pack
No.5 (24 mm)
7312 Pastel Yellow
7313 Pastel Blue
7314 Pastel Mint
7315 Pastel Lavender ...
7317 Pastel Pink

HAHNENKRATT BLACKLINE ENDO MIRRORS ULTRA FS TITANIUM FRONT SURFACE BENDABLE AUTOCLAVABLE



- Micro mirrors with bendable stem for a precise alignment of the mirror
- Can be also used for Periodontics especially for flap openings

BLACKline Microflex Endo Mirrors
Titanium Front Surface Bendable
B5181 ULTRA FS Round 3 mm
B5182 ULTRA FS Round 5 mm
B5183 ULTRA FS Oval 3 x 6 mm
B5184 ULTRA FS Drop Shaped 3 x 6 mm
B442 Ergotouch Black Mirror Handle Stem ..

HAHNENKRATT BLACKLINE MIRRORS ULTRA FS TITANIUM FRONT SURFACE AUTOCLAVABLE





- DLC coating (Diamond-Like Carbon) for a high-end finish
- Anti reflective, dark black coating
- Avoids visual fatigue by working under light and magnification


BLACKline Front Surface Titanium ULTRA FS Mirrors Stem
B6160 ULTRA FS Single Side, 3 pcs/pack No.4 (22 mm)
B6161 ULTRA FS Single Side, 3 pcs/pack No.5 (24 mm)
B6640 ULTRA FS Double Sides, 3 pcs/pack No.4 (22 mm) .
B6641 ULTRA FS Double Sides, 3 pcs/pack No.5 (24 mm) .
B442 Ergotouch Black Mirror Handle Stem




AMALGAM INSTRUMENTS HOLLOW HANDLE EXTRA LIGHT



FILLING INSTRUMENTS HOLLOW HANDLE EXTRA LIGHT



EXCAVATORS HOLLOW HANDLE EXTRA LIGHT



AUTOCLAVABLE IMPLANT SCALER SET





- Fully autoclavable up to 134°C
- Replaceable carbon reinforced resin tips, can be resharpened similar to metal tips
- Extra light hollow handle

1128 Implant Scaler H6/H7 Hollow Handle
Handle + 24 Tips (12 each H6/H7) ...
24 Tips (12 each H6/H7)

MIRROR HANDLES HOLLOW EXTRA LIGHT



SCALERS AND CURETTES HOLLOW HANDLE EXTRA LIGHT




PROBES HOLLOW HANDLE EXTRA LIGHT




PERIO PROBES LASER MARKED HOLLOW HANDLE EXTRA LIGHT

DOUBLE ENDED



SINGLE ENDED




EXAMINATION INSTRUMENTS STAINLESS STEEL




EXCAVATORS STAINLESS STEEL



FILLING INSTRUMENTS STAINLESS STEEL



AMALGAM INSTRUMENTS STAINLESS STEEL



CEMENT SPATULAS STAINLESS STEEL



AMALGAM WELL STAINLESS STEEL



53-030-001
Amalgam well


RUBBER DAM PUNCHES AND FORCEPS



EDENTA MANDRELS STAINLESS STEEL



EXTRACTION FORCEPS ANATOMICAL HANDLE



SURGICAL INSTRUMENTS STAINLESS STEEL



SPLINTER FORCEPS STAINLESS STEEL


TWEEZERS STAINLESS STEEL




PROBES STAINLESS STEEL



ROOT ELEVATORS LUXATING TYPE STAINLESS STEEL



ROOT ELEVATORS STAINLESS STEEL



SCISSORS

STAINLESS STEEL



62-210-115
Iris Straight
Smooth 11 cm



62-211-115
Iris Curved
Smooth 11 cm

TUNGSTEN CARBIDE



10-084-115
Iris Straight
Smooth 11.5 cm



10-085-115
Iris Curved
Smooth 11.5 cm



10-360-140
Hegar
14 cm
€3.9



10-414-140
Mathieu
14 cm
€5.9

NEEDLE HOLDERS TUNGSTEN CARBIDE TIPS


NEEDLE HOLDERS STAINLESS STEEL



69-124-135
Hegar
13.5 cm



69-124-150
Hegar
15 cm



69-213-140
Mathieu
14 cm

HEMOSTATIC FORCEPS STAINLESS STEEL



66-054-125
Mosquito
Straight 12 cm



66-055-125
Mosquito
Curved 12 cm




66-070-140
Kelly
Straight 14 cm



66-071-140
Kelly
Curved 14 cm

SEMKEN TISSUE FORCEP



64-191-120
Tissue Forcep
12 cm



actual size

INJECTION SYRINGES



60-820-001 Standard 2.2 ml ...



60-855-001 Aspirating 2.2 ml ...

CHEEK RETRACTORS STAINLESS STEEL



20-121-140 Sternberg 14 cm
20-121-160 Sternberg 16 cm
20-122-145 Roux-Magnus 14.5 cm ..

CURVED MIXING SPATULA AUTOCLAVABLE



- Light weight 58 gm only
- Fully autoclavable up to 134°C
- Stainless steel blade
- Handle colour: Blue

455
Curved Spatula

ALGinate AND PLASTER SPATULAS AUTOCLAVABLE



450XS
Extra Stiff Grey

450R
Regular Blue

- Fully autoclavable up to 134°C
- Fit all flexible mixing bowls
- Flexible control in mixing with easy grip handle
- Will not scratch the surface of bowls and clean up easily

MATRIX RETAINERS STAINLESS STEEL



58-030-005 Ivory ...



58-020-001 Universal ...



58-020-000 Junior ...

MICRODONT TOFFLEMIRE MATRIX BANDS STAINLESS STEEL



10.405.001
Tofflemire Matrix Band No. 1
12 pcs/pack



10.405.002
Tofflemire Matrix Band No. 2
12 pcs/pack

PLASTIC AMALGAM CARRIERS AUTOCLAVABLE



A11 Plastic Amalgam Gun Curved ...



A12 Plastic Amalgam Gun Straight ...

- Autoclavable up to 134°C
- Replacement tips A11 ...
- Replacement tips A12 ...



IMPRESSION TRAYS STAINLESS STEEL PERF-LOCK

UPPER	NO.	WIDTH (mm)	DEPTH (mm)
Extra-Large	60-120-001	76	64
Large	60-120-002	71	60
Medium	60-120-003	71	58
Small	60-120-004	65	54
Extra-Small	60-120-005	65	52



LOWER	NO.	WIDTH (mm)	DEPTH (mm)
Extra-Large	60-130-001	75	57
Large	60-130-002	75	54
Medium	60-130-003	72	54
Small	60-130-004	62	50
Extra-Small	60-130-005	55	45



PERF-PARTIAL	NO.	WIDTH (mm)	DEPTH (mm)
Upper Anterior	61-180-004	66	38
Lower Anterior	61-180-003	64	36
Quadrant	61-180-002	32	54
Quadrant	61-180-001	32	54

60-56-123 Adult Set, 8 pcs/set ...
Individual Tray



IMPRESSION TRAYS STAINLESS STEEL PERF-LOCK EDENTULOUS

UPPER	No.	WIDTH (mm)	DEPTH (mm)
Extra-Large	60-160-001	76	70
Large	60-160-002	73	66
Medium	60-160-003	70	66
Small	60-160-004	68	64

LOWER	No.	WIDTH (mm)	DEPTH (mm)
Extra-Large	61-170-001	82	62
Large	61-170-002	70	62
Medium	61-170-003	68	60
Small	61-170-004	70	56

60-56-456 Edentulous Set, 8 pcs/box ...
Individual Tray

MIXING SPATULAS AND PLASTER KNIVES



WAXING UP INSTRUMENTS STAINLESS



SUNDRY INSTRUMENTS STAINLESS STEEL



TUNGSTEN CARBIDE TIPS PLIERS



€6.85



HAHNENKRATT FRONT SURFACE MIRRORS ULTRA FS AND ULTRAvision FS

The split mirror shows
the crystal clear bright
reflection of the **ULTRA FS**
compared to Rhodium

ULTRA FS and **ULTRAvision FS**:- The brightest and sharpest mirrors
worldwide. The new unique design of the mirror layers provides
superior bright images and crystal clarity previously unavailable

- With about 113% reflection factor, these
ULTRA FS and **ULTRAvision FS** mirrors
are even brighter than Rhodium
- Superior colour accuracy
- Exact sharpness
- Scratch resistant

FRONT SURFACE TITANIUM MIRRORS DOUBLE SIDED



- Available in sizes 4 and 5
- Helps illuminate, retract and control tongue
- Autoclavable to **134°C**
- Front surface reflecting
- Simple stem fitting

7624
Titanium Double Sided Mirrors
No.5 **24 mm**
6 pcs/box

HAHNENKRATT FRONT SURFACE MIRRORS MEGA RHODIUM



- Frame and stem made out of one single
piece of stainless steel, considering
ergonomic aspects; - this characterizes
the uniqueness of the MEGA
- Form follows function - it's ergonomically
designed

MEGA Rhodium Stem
6060, No. 4 **22 mm**, 12 pcs/box ..
6061, No. 5 **24 mm**, 12 pcs/box ..

MEGA DOUBLE SIDED Rhodium Stem
6100, No. 4 **22 mm**, 6 pcs/box
6101, No. 5 **24 mm**, 6 pcs/box

HANDLE NOT INCLUDED



No.3 (20 mm)
No.4 (22 mm)
No.5 (24 mm)

MOUTH MIRRORS STAINLESS STEEL

Mouth Mirrors Stem, 10 pcs/box
47060020 Plain No.3
47060022 Plain No.4
47060024 Plain No.5
47062020 Magnifying No.3
47062022 Magnifying No.4
47062022 Magnifying No.5
47065020 Front Surface No.3 ...
47065022 Front Surface No.4 ...
47065024 Front Surface No.5 ...

MIRROR HANDLES STAINLESS STEEL



47-050-003 Mirror Handle Round Stem ...



47-050-004 Mirror Handle Octagonal Stem ...



47-050-008 Mirror Handle with Rule ...



47-051-002 Hollow Mirror Handle **7.5 mm** Extra Light ...

MIRROR HANDLES HOLLOW EXTRA LIGHT



HH-80-006 Hollow Mirror Handle **XL Stem** ..



HH-80-005 Hollow Mirror Handle **XL Titanium Coated Stem** ...

New Edition

